

Product description:

The software is designed for creating formulas for measured mixtures with possibility of declaring the tolerance field for particular ingredient and measuring each ingredient in the receipt.



Nowa		SCALER: 1.575	
0,0 CAŁKOWITA MASA SKŁADNIKA		NAZWA RECEPTURY: Nowa	
000,000 ODCZYT Z WAGI		IŁOŚĆ SZARŻ: 10	DOPUSZCZALNE ODCHYLENIE: 0%
-100,0% ODCHYLENIE RZECZYWISTE		NR SZARŻY: 1	
2,00025 NORMA		1 INGREDIENT: BIOTYNA (H)	
TOO FEV		KOD: _____	
PANEL STEROWANIA			
NASTĘPNY SKŁADNIK		WYJŚCIE	
ANULUJ		ZMIANA FORMALKI	

Designed for scales:

- all weighing instruments offered by RADWAG (in single stand version)

Functions:

- database for formulas and ingredients located in the software
- executing the weighing processes of recipes
- creating reports from accomplished measurements
- measuring a formula for a particular ingredient
- measuring a formula to determined mass
- reporting of an indication of material consumption

Minimal hardware requirements:

- Processor 1,2 GHz
- Minimum of HDD 500 MB free space
- RAM 256 MB (reccommended 512 MB)
- OS Windows 2000/XP

* NOTE! RADWAG does not sell software "in box". Boxes presented on our Web-site are of virtual characteristic only, and they are not accessible with purchased software. The customer receives purchased software on a CD in special protective cover and the user manual as a brochure.