



Foot button is designed to send to a weighing instruments commands on TARRING or PRINTING (print measurement result).

Application of the foot tare and print button enables carrying out complex procedures of weighing, checkweighing, filling in cases, where both operator's hands are involved in the procedure execution. Pressing foot button allows for operating the balance with hands and simultaneous tarring, printing and sending the result to a computer. The cooperation of the system comprising a weighing instrument - computer - PRINT foot button - TARE foot button is ensured by signal distributor which is part of the foot buttons set.

Set of TARE or PRINT buttons is applicable for laboratory weighing instrument by RADWAG equipped with graphic display and RADWAG industrial scales featuring indicators or terminals series: TMC and PUE C41/H.

