


RADWAG Laboratory Holders



Discover new RADWAG Laboratory Holders designed for faster and more comfortable weighing!

This is perfect and secure fixing of small vessels both typical and non-standard shapes.





Using RADWAG Laboratory Holders ensures full ergonomics of weighing.

Laboratory holders provide dosing option of a sample, which is often of high value, directly into a vessel.

Thanks to that the efficiency is increased, costs are reduced and loss resulting from moving samples is minimized.

Holders have been designed to facilitate dosing samples directly into typical vessels such as: measuring flasks, round flasks, test tubes, titration vessels.

Laboratory holders can be used with all Radwag's analytical balances.

UL-01	UL-02	UL-03	UL-04
 <p>Purpose:</p> <ul style="list-style-type: none"> • test tubes • filters • laboratory glassware 	 <p>Purpose:</p> <ul style="list-style-type: none"> • test tubes ø 8,10,12 mm • centrifuge tubes ø 8,10,12 mm 	 <p>Purpose:</p> <ul style="list-style-type: none"> • round flasks of 50, 100, 250 ml 	 <p>Purpose:</p> <ul style="list-style-type: none"> • flat flasks of 50, 100, 250 ml
For: AS/2Y, AS/X, AS	For: AS/2Y, AS/X, AS	For: AS/2Y, AS/X, AS	For: AS/2Y, AS/X, AS
